

Common Questions Concerning the Pharmanex® BioPhotonic Scanner

Question:

What science is there to validate the Pharmanex BioPhotonic Scanner?

RESPONSE:

The Pharmanex BioPhotonic Scanner is backed by science.

The use of Raman spectroscopy for biological measurements is an established scientific discipline backed by years of research. The Pharmanex BioPhotonic Scanner is a patented application of Raman spectroscopy for the measurement of carotenoid antioxidant nutrients in living tissue for the improvement of nutrition. The use of biophotonics to assess biological molecules in living tissue is a distinct scientific discipline, and the Pharmanex BioPhotonic Scanner is an instrument that is based on this scientific discipline.

The use of Raman spectroscopy for the assessment of human tissue carotenoids has been validated by at least eight peer-reviewed studies conducted by third party entities unrelated to Pharmanex or the supplementation industry. (Bernstein, 1998, 2002; Ermakov, 2004a, 2004b; Gellermann, 2004, 2002; Hata, 2000; Zhao, 2003). Raman spectroscopy is an established and accepted detection method; a Medline search for the term *raman spectroscopy* yields over 4,400 articles. Most importantly, validation of Raman spectroscopy is recognized through the award of Nobel Prize to Sir CV Raman for its discovery in 1930.

Pharmanex is not the only entity impressed with the nutritional implications of Raman spectroscopy. In 2003 the National Cancer Institute (a division of the National Institutes of Health) awarded researchers at Yale University a \$1 million grant to conduct a study using Raman detection of carotenoids as an objective measure of fruit and vegetable intake. Most large-scale nutritional studies rely on diet recall surveys, which are subject to errors of reporting fruit and vegetable consumption. The study will be completed in 2006 (Mayne *et al*, 2003; YCC News Release, 2003).

Importantly, Pharmanex has validated the use of Raman spectroscopy for the measurement of carotenoids in 3 studies including a large-scale clinical screening study with 1,375 subjects that confirmed a correlation between antioxidant status and lifestyle parameters (Smidt, 2003). A second study established efficacy of LifePak® to improve the antioxidant status of subjects over a 12-week period (Smidt, 2002), and a third study established a highly significant correlation ($r=0.78$, $p < 0.001$) between serum carotenoid levels and skin carotenoid levels as assessed by the Pharmanex BioPhotonic scanner (Smidt, 2004a). More studies are ongoing; a study in conjunction with the University of Utah looking at the correlation between skin carotenoid levels and levels of water soluble antioxidants; a large-scale Pharmanex in-house correlation study; and a study conducted at the Rippe Lifestyle Institute focusing on the antioxidant effects of weight-loss. These studies after peer review will be published in scientific journals.

The Pharmanex BioPhotonic Scanner is consistently well received by experts in all areas of science. Pharmanex Scientists have presented the science behind LifePak® and the BioPhotonic scanner at a number of scientific meetings. In February 2003, Pharmanex scientists joined with Dr. James Rippe and the inventor of the biophotonic scanner, Dr. Werner Gellermann in a presentation at the New York Academy of Sciences. In February 2003 and 2004, Pharmanex Scientists attended scientific meetings on antioxidants organized by the Oxygen Club of California in Cadiz, Spain, and Santa Barbara, California. Dr. Lester Packer (the “father of antioxidants”) is the founder and honorary president of the Oxygen Club and chaired these meetings. Dr. Carsten Smidt presented a scanner study at the Santa Barbara meeting. In April 2003 and 2004, Dr. Carsten Smidt attended the Federation of American Societies for Experimental Biology (FASEB) meetings in San Diego and Washington DC, and presented the results of two different scanner studies. The FASEB meeting is attended by more than 10,000 scientists from around the world, but only a small percentage are selected to present. In January of 2004, Dr. Carsten Smidt attended and presented the scanner/serum correlation study at the prestigious Gordon Research Conference on Carotenoids in Ventura, California. Dr. Smidt also scanned 60 of the top antioxidant researchers in the world at this conference. In all instances the Pharmanex BioPhotonic Scanner has been very well received by the scientific community.

Question:

If the BioPhotonic Scanner measures only carotenoids, how can it be used to infer the status of other nutrients?

RESPONSE:**Serum carotenoids correlate to overall antioxidant status**

Carotenoid molecules are not regenerated like other antioxidants, and are degraded in the process of neutralizing free radicals or reactive oxygen species. A typical carotenoid molecule like lycopene or β -carotene is able to sustain more than 20 free radical hits by lipid radicals before it becomes completely destroyed (Tsuchiya, 1994). Lycopene and β -carotene are just two examples of antioxidants among hundreds of antioxidants that make up the *antioxidant network*. Carotenoids act sacrificially to protect other members of the antioxidant network (such as vitamins E and C) from having to sustain free radical hits; in this way carotenoids will support the entire antioxidant network consequently reducing the danger from oxidative stress (Packer, 1994; Packer and Coleman, 1999). Conversely, high levels of oxidative stress (e.g., with smoking) adversely affect the antioxidant network, and the resulting increased free radical activity leads to a depletion or reduction in tissue carotenoids (Smidt and Shieh, 2003; Gollnick and Siebenwirth, 2002; Dietrich, 2003).

A recent study conducted by Svilaas *et al.* established carotenoids as a reliable indicator of other dietary antioxidants. Svilaas and his colleagues assessed antioxidant intake from diets of more than 2,670 adults, and evaluated blood serum antioxidants of 61 individuals for seven consecutive days. Svilaas *et al.* found the ability of carotenoids to predict serum levels of other antioxidants was stronger than the predictive ability of alpha, beta, delta, and gamma-tocopherols as well as glutathione (Svilaas, 2004). Carotenoids are not only convenient biomarkers because they are accurate predictors of overall antioxidant status, but also they are Raman active and can be detected without the concerns of blood samples (Bernstein, 1998, 2002; Gellermann 2002a; Zhao, 2003; Ermakov, 2004b). Furthermore, carotenoids are delivered to tissues by the same mechanism as other fat-soluble antioxidants. This shared LDL delivery mode is the proposed mechanism to explain the correlation between tissue carotenoids and other fat-soluble antioxidants in multiple studies (Lasheras *et al.*, 2002; Steinberg & Chait, 1998).

Serum carotenoids correlate to skin carotenoids

Skin carotenoids analyzed by HPLC were shown to correlate significantly to serum carotenoid levels (Peng, 1995). The significant correlation between skin biopsy-levels of carotenoids and serum carotenoid levels eliminates the need for routine skin removal. A more recent study showed a highly significant correlation between serum total carotenoids and skin carotenoids as assessed by Raman Spectroscopy ($r = 0.78$, $p < 0.001$) (Smidt, 2004). These data bridge the findings of Svilaas, and Peng, to validate Raman Spectroscopy as a method to assess skin carotenoid status as an indication of broad-spectrum antioxidant status, without the inconvenience of skin and blood samples.

Carotenoids accepted as indicator of fruit and vegetable intake

Carotenoids include more than 50 antioxidants widely distributed among fruits and vegetables. When ingested from dietary sources, the presence of carotenoids in living tissue is an indicator for the presence of other important nutrients common to those dietary sources (Svilaas, 2004). Based on this correlation, the National Cancer Institute awarded a \$1 million research-grant to Yale scientists to conduct a study using the Raman-detection of carotenoids as an objective measure of fruit and vegetable intake (Mayne *et al.*, 2003; YCC News Release, 2003). The correlation between carotenoids and other nutrients applies to nutritional supplements only to the degree they deliver optimal amounts of all essential and generally beneficial nutrients. For this reason, LifePak[®] is formulated as a broad-spectrum multivitamin, mineral, antioxidant supplement, which contains nutrients in amounts similar to a diet rich in fruits and vegetables.

Question:

What aspects of health are positively affected by proper carotenoid nutrition?

RESPONSE:

Carotenoids has been shown in countless studies to support many areas of health

The scanner is not intended to diagnose, mitigate, treat, or cure any disease. Nonetheless, convincing evidence suggests that certain carotenoids have been linked to health benefits including reduced risk of age-related macular degeneration, cataracts, cardiovascular disease, and prostate cancer. A review article written by Pharmanex scientists will appear in the peer-reviewed journal *Current Trends in Nutraceutical Research*. The article includes a review of the role of carotenoids in human health and is summarized below (Smidt and Burke, 2004).

Eye Health

A number of studies support the protective role of carotenoids in the prevention of age-related eye diseases. For example, reduced risks of cataracts (Brown *et al.*, 1999; Chasan-Taber *et al.*, 1999) and age-related macular degeneration (Seddon, 1994; Bone, 2000, 2003; Landrum, 1997, 1996; Elless, 2000; Bernstein, 2002, Richer, 1999; Hammond, 1997) have been associated with high intakes of vegetables rich in the carotenoids lutein and zeaxanthin.

Cardiovascular Health

The carotenoids lutein and lycopene have been shown separately to support multiple aspects of cardiovascular health. Studies show positive health implications of lutein including decreased risk of mortality from cardiovascular disease (Kouris-Blazos, 2002), decreased progression of pre-atherosclerotic conditions (Dwyer, 2001), and other cardioprotective effects (Olmedilla, 2001; Cardinault, 2003).

The cardioprotective effects of lycopene have also been shown in multiple studies including reduced risk of myocardial infarction (Kohlmeier, 1997), lower risk of cardiovascular disease (Sesso, 2004), reduced LDL oxidation (Agarwal, 1998), and reduced production of LDL cholesterol (Fuhrman *et al.*, 1997).

Cancer

Epidemiological studies have shown that high intakes of tomatoes and tomato products, rich in lycopene, as well as high blood levels of lycopene are significantly associated with decreased prostate cancer risk (Deming, 2002; Giles, 1997; Giovanucci, 1995, 2002; Lu, 2001; Vogt, 2002). The finding of Kucuk *et al.* suggests that lycopene supplementation may decrease prostate cancer growth (Kucuk *et al.*, 2001). These effects may be attributed to lycopene's antioxidant and DNA protective properties (Riso *et al.*, 1999; Porrini and Riso, 2000).

Carotenoids may also play a role in cancer prevention because they can enhance gap junctional communication (GJC) between cells. (Krutovskikh *et al.*, 1995; Yamasaki *et al.*, 1995; Yamasaki, 1995; Dahl *et al.*, 1995; Trosko, 2003). Lycopene and β -carotene have been shown to enhance GJC significantly, and these effects are not related to their known antioxidant properties (Sies and Stahl, 1997; Stahl *et al.*, 1997; Zhang *et al.*, 1991). Thus, carotenoids may act via two distinct mechanisms of action to protect from cancer: as antioxidants to prevent mutagenic DNA alterations, and as promoters of GJC.

Oxidative Stress correlates to skin carotenoid concentration

Many of the above diseases are known to be linked to oxidative stress. A population study of 1,375 subjects was conducted at the Pharmanex Research Institute and found that individuals with high oxidative stress generally have low skin carotenoid levels as measured by Raman spectroscopy, independent of subjects' dietary carotenoid consumption. This correlation was demonstrated by using Urinary MDA test, a proven model for oxidative stress (Smidt and Shieh, 2003).

Question:

Are other methods of testing more accurate at assessing antioxidant status?

RESPONSE:

Raman spectroscopy has been shown to be more reliable than other methods and more suitable for routine measurements reduced probability of error

Before Raman detection of carotenoids, established methods of measuring antioxidant status included the analysis of blood, urine, or tissue samples. These tests are expensive, invasive (e.g. *skin biopsies; blood draw*), have a higher probability of error (multiple steps involved in sample preparation and quantification all increase the magnitude of the overall error). Moreover, it may take weeks to receive the results. Also,

the accuracy of tests using blood/serum or urine is placed in question due to possible effects of recent meals as the timing that the biological sample is taken (blood or urine) influences the outcome. A recent clinical study showed that Raman measurement of carotenoids in humans skin (palm of hand) correlates very highly with total serum carotenoid levels quantified from fasted serum (Smidt, Gellermann, & Zidichouski, 2004). The BioPhotonic Scanner is a great way to measure carotenoid levels safely and non-invasively at the site of action with the added advantage that it is well suited for routine measurements in large populations. The ability to assess carotenoid status through Raman spectroscopy will lead to further advances in nutritional and biological measurement.

Question:

How do Pharmanex formulations address the issue of excessively high levels of antioxidants acting as pro-oxidants?

RESPONSE:

Pharmanex products contain nutrient amounts proven to be safe in clinical studies

The levels of all nutrients found in Pharmanex products are based on well documented epidemiological, clinical, pre-clinical and safety studies. LifePak[®] is formulated to provide optimal nutrition with substantiated levels of nutrients that will not induce a pro-oxidative state. Included in the comprehensive blend of antioxidants is a balanced carotenoid combination in amounts similar to those provided by diets high in fruits and vegetables: 7.5 mg β -carotene, 5 mg lycopene, 2 mg α -carotene and 2 mg lutein. Each ingredient in LifePak[®] is present in amounts that are documented to be safe for long-term supplementation. Further, LifePak is safe when taken in conjunction with a diet high in fruits and vegetables. The daily amounts of all vitamins and minerals are well below the No-Observed Adverse Effect Levels (NOAEL) established by the Council for Responsible Nutrition (CRN) in 1997 and the Upper Limits (UL) established by the Food and Nutrition Board of the National Research Council.

BIBLIOGRAPHY

Agarwal ,S. and Rao ,A.V. (1998) Food and nutrient intakes of individuals in 1 day in the United States. Preliminary Report #2. 1980. Tomato lycopene and low density lipoprotein oxidation: A human dietary intervention study. *Lipids* **33**, 981-984.

Bernstein P.S., Raman detection of macular carotenoid pigments in intact human retina. *Invest Ophthalmol Vis Sci.* 1998 Oct;39(11):2003-11.
http://www.ncbi.nlm.nih.gov/entrez/query.fcgi?cmd=Retrieve&db=pubmed&dopt=Abstract&list_uids=9761278

Bernstein, P.S. et al, Resonance Raman measurement of macular carotenoids in normal subjects and in age-related macular degeneration patients. *Ophthalmology.* 2002 Oct;109(10):1780-7.
http://www.ncbi.nlm.nih.gov/entrez/query.fcgi?cmd=Retrieve&db=pubmed&dopt=Abstract&list_uids=12359594

Bernstein ,P.S., Gellermann W. Measurement of carotenoids in the living primate eye using resonance Raman spectroscopy. *Methods Mol Biol.* 2002;196:321-9.
<http://biomed.humanapress.com/ChapterDetail.pasp?isbn=1-59259-274-0&ccode=1-59259-274-0:321>

Bone, R.A., Landrum, J.T., Dixon, Z., Chen, Y. and Llerena,C.M. (2000) Lutein and zeaxanthin in the eyes, serum and diet of human subjects. *Experimental Eye Research* **71**, 239-245.

Bone, R.A., Landrum, J.T., Mayne, S.T., Gomez, C.M., Tibor, S.E. and Twaroska, E.E. (2001) Macular pigment in donor eyes with and without AMD: a case-control study. *Investigative Ophthalmology and Visual Science* **42**, 235-240.

Brown, L., Rimm, E.B., Seddon, J.M., Giovanucci, E.L., Chasan-Taber, L., Spiegelman, D., Willett, W.C. and Hankinson, S.E. (1999) A prospective study of carotenoid intake and risk of cataract extraction in US men. *American Journal of Clinical Nutrition* **70**, 517-524.

Cardinault, N., Gorrard, J.M., Tyssandier, V., Grolier, P., Rock, E. and Borel, P. (2003) Short-term supplementation with lutein affects biomarkers of lutein status similarly in young and elderly subjects. *Experimental Gerontology* **38**, 573-582.

Chasan-Taber, L., Willett, W.C., Seddon, J.M., Stampfer, M.J., Rosner, B., Colditz, G.A., Speizer, F.E. and Hankinson, S.E. (1999) A prospective study of carotenoid and vitamin A intakes and risk of cataract extraction in US women [see comments]. *American Journal of Clinical Nutrition* **70**, 509-516.

Dahl, E., Winterhager, E., Traub, O. and Willeck, K (1995) Expression of gap junction genes, connexin40 and connexin43, during fetal mouse development. *Anatomy and Embryology* **191**, 267-278.

Dietrich M.,Block G.,Hudes M.,Morrow J.D.,Norkus E.P.,Traber M.G.,Cross C.E. and Packer L. Antioxidant supplementation decreases lipid peroxidation biomarker F (2)-isoprostanes in plasma of smokers. *Cancer Epidemio Biomarkers Prev* 2002;1:7-13.

Deming, D.M., Boileau, T.W-M, Heintz, K.H., Atkinson, C.A. and Erdman, J.W., Jr. (2002) Carotenoids: Linking chemistry, absorption, and metabolism to potential roles in human health and disease. In: Cadenas, E. and Packer, L. (Eds), *Handbook of Antioxidants* (New York: New York: Marcel-Dekker.), pp. 189-221.

Dwyer, J.H., Navab, M., Dwyer, K.M., Hassan, K., Sun, P., Shircore, A., Hama-Lavy, S., Hough, G., Wang, X., Drake, T., Merz, C.N. and Fogelman, A.M. (2001) Oxygenated carotenoid lutein and progression of early atherosclerosis: the Los Angeles atherosclerosis study. *Circulation* **103**, 2922-2927.

Elles, M.P., Blaylock, M.J., Huang, J.W. and Gussman, C.D. (2000) Plants as a natural source of concentrated mineral nutritional supplements. *Food Chemistry* **71**, 181-188.

Ermakov, I.V. et al. Noninvasive selective detection of lycopene and beta-carotene in human skin using Raman spectroscopy. *J Biomed Opt.* 2004a Mar;9(2):332-8.
http://www.ncbi.nlm.nih.gov/entrez/query.fcgi?cmd=Retrieve&db=pubmed&dopt=Abstract&list_uids=15065899

Ermakov I.V. et al. Macular pigment Raman detector for clinical applications. *J Biomed Opt.* 2004b Jan-Feb;9(1):139-48.
http://www.ncbi.nlm.nih.gov/entrez/query.fcgi?cmd=Retrieve&db=pubmed&dopt=Abstract&list_uids=14715066

Fuhrman, B., Elis, A. and Aviram, M. (1997) Hypocholesterolemic effect of lycopene and beta-carotene is related to suppression of cholesterol synthesis and augmentation of LDL receptor activity in macrophages. *Biochemical and Biophysical Research Communications* **233**, 658-662.

Gellermann, W., Bernstein PS. Noninvasive detection of macular pigments in the human eye. *J Biomed Opt.* 2004 Jan-Feb;9(1):75-85. Review.
http://www.ncbi.nlm.nih.gov/entrez/query.fcgi?cmd=Retrieve&db=pubmed&dopt=Abstract&list_uids=14715058

Gellermann W. et al, In vivo resonant Raman measurement of macular carotenoid pigments in the young and the aging human retina. *J Opt Soc Am A Opt Image Sci Vis.* 2002a Jun;19(6):1172-86.
http://www.ncbi.nlm.nih.gov/entrez/query.fcgi?cmd=Retrieve&db=pubmed&dopt=Abstract&list_uids=12049355

Giles, G., Ireland, P. (1997) Diet, nutrition and prostate cancer. *International Journal of Cancer* **72**, 13-17.

Giovanucci, E., Ascherio, A., Rimm, E.B., Stampfer, M.J., Colditz, G.A. and Willett, W. (1995) Intake of carotenoids and retinol in relation to risk of prostate cancer. *Journal of the National Cancer Institute* **87**, 1767-1776.

Giovanucci, E., Rimm, E.B., Stampfer, M.J. and Willett, W.C. (2002) A prospective study of tomato products, lycopene and prostate cancer risk. *Journal of the National Cancer Institute* **94**, 391-398.

Gollnick HP, Siebenwirth C. Beta-carotene plasma levels and content in oral mucosal epithelium is skin type associated. *Skin Pharmacol Appl Skin Physiol*. 2002 Sep-Oct;15(5):360-6.

Hammond, B.R., Jr., Johnson, E.J., Russell, R.M., Krinsky, N.I., Yeum, K.J., Edwards, R.B. and Snodderly, D.M. (1997) Dietary modification of human macular pigment density. *Investigative Ophthalmology and Visual Science* **38**, 1795-1801.

Hata, T.R. et al, Non-invasive raman spectroscopic detection of carotenoids in human skin. *J Invest Dermatol*. 2000 Sep;115(3):441-8.
http://www.ncbi.nlm.nih.gov/entrez/query.fcgi?cmd=Retrieve&db=pubmed&dopt=Abstract&list_uids=10951281

Kohlmeier, L., Kark, J.D., Gomez-Gracia, E., Martin, B.C., Steck, S.E., Kardinaal, A.F., Ringstad, J., Masaev, V., Riemersma, R., Martin-Moreno, J.M., Huttunen, J.K. and Kok, F.J. (1997) Lycopene and myocardial infarction risk in the EURAMIC Study. *American Journal of Epidemiology* **146**, 618-626.

Kouris-Blazos, A. (2002) Morbidity mortality paradox of 1st generation Greek Australians. *Asia Pacific Journal of Clinical Nutrition* **11 Suppl 3**, S569-S575.

Kucuk, O., Sarkar, F.H., Sark, W., Djuric, Z., Pollak, M.N., Khachik, F., Li, Y.W., Banerjee, M., Grignon, D., Bertram, J.S., Crissman, J.D., Pontes, E.J. and Wood, D.P., Jr. (2001) Phase II randomized clinical trial of lycopene supplementation before radical prostatectomy. *Cancer Epidemiology Biomarkers & Prevention* **10**, 861-868.

Krutovskikh, V.A., Mesnil, M., Mazzoleni, G. and Yamasaki, H. (1995) Inhibition of rat liver gap junction intercellular communication by tumor-promoting agents in vivo. Association with aberrant localization of connexin proteins. *Laboratory Investigations* **72**, 571-577.

Lasheras C, Huerta JM, Gonzalez S, Brana AF, Patterson AM, Fernandez S. Independent and interactive association of blood antioxidants and oxidative damage in elderly people. *Free Radic Res*. 2002 Aug;36(8):875-82.

Landrum, J.T., Bone, R.A., Joa, H., Kilburn, M.D., Moore, L.L. and Sprague, K.E.. (1997) A one year study of the macular pigment: The effect of 140 days of a lutein supplement. *Experimental Eye Research* **65**, 57-62.

Landrum, J.T., Bone, R.A., Kilburn, M.D., Joa, H. and Gomez, C (1996) Dietary lutein supplementation increases macular pigment (MP). *Faseb Journal* **10**:A242 (Abstract).

Lu, Q.Y., Hung, J.C., Heber, D., Liang, V., Go, W., Reuter, V.E., Cordon-Cardo, D., Scher, H.I., Marshall, J.R. and Zhang, Z.F. (2001) Inverse associations between plasma lycopene and other carotenoids and prostate cancer. *Cancer Epidemiology Biomarkers & Prevention* **10**, 749-756.

Mayne, S.T., NIH funded study in progress: Novel, Noninvasive Biomarker of Fruit & Vegetable Intake, Computer Retrieval of Information on Scientific Projects. Grant Number 1R01CA096838-01A1.
http://crisp.cit.nih.gov/crisp/CRISP_LIB.getdoc?textkey=6610240&p_grant_num=1R01CA096838-01A1&p_query=&ticket=8698486&p_audit_session_id=39073479&p_keywords=

Olmedilla, B., Granado, F., Southon, S., Wright, A.J., Blanco, I., Gil-Martinez, E., Berg, H., Corridan, B., Roussel, A.M., Chopram, M. and Thurnham, D.I. (2001) Serum concentrations of carotenoids and vitamins A, E, and C in control subjects from five European countries. *British Journal of Nutrition* **85**, 227-238.

Packer L. Vitamin E is nature's master antioxidant. *Scientific American, Science and Medicine* 1994, 1:54-63.

Packer L, Coleman C. *The Antioxidant Miracle*. John Wiley and Sons, New York, 1999.

Peng, Y.M., Peng, Y.S., Lin, Y., Moon, T., Roe, D.J. and Ritenbaugh, C. (1995) Concentrations and plasma-tissue-diet relationships of carotenoids, retinoids, and tocopherols in humans. *Nutrition and Cancer* **23**, 233-246.

Porrini, M. and Riso P. (2000) Lymphocyte lycopene concentration and DNA protection from oxidative damage is increased in women after a short period of tomato consumption. *Journal of Nutrition* **130**, 189-192.

Richer, S. (1999) ARMD--pilot (case series) environmental intervention data [In Process Citation]. *Journal of the American Optometric Association* **70**, 24-36.

Riso, P., Pinder, A., Santangelo, A. and Porrini, M. (1999) Does tomato consumption effectively increase the resistance of lymphocyte DNA to oxidative damage? *American Journal of Clinical Nutrition* **69**, 712-718.

Sies, H. and Stahl, W. (1997) Carotenoids and intercellular communication via gap junctions. *International Journal for Vitamin and Nutrition Research* **67**, 364-367.

Seddon, J.M., Ajani, U.A., Sperduto, R.D., Hiller, R., Blair, N., Burton, T.C., Farber, M.D., Gragoudas, E.S., Haller, J., Miller, D.T. Yannuzzi, L.A. and Willett, W. (1994) Dietary carotenoids, vitamins A, C, and E, and advanced age-related macular degeneration. Eye Disease Case-Control Study Group. *Journal of the American Medical Association* **272**, 1413-1420.

Sesso, H.D., Buring, J.E., Norkus, E.P., and Gaziano, J.M. (2004) Plasma lycopene, other carotenoids and retinol and the risk of cardiovascular disease in women. *American Journal of Clinical Nutrition* **79**, 47-53.

Smidt, C.R., Gellermann, W., Zidichouski, J.A. Non-invasive Raman spectroscopy measurement of human carotenoid status. Pharmanex Research Institute. FASEB 18(4):A480, 2004a.

Smidt, C.R., Burke, D.S. Nutritional Significance and Measurement of Carotenoids. *Curr Trends Nutr Res*. To be published in May 2004b issue.

Smidt, C.R., Clinical Screening Study: Use of the Pharmanex BioPhotonic Scanner to assess skin carotenoids as a marker of antioxidant status. Pharmanex in-house Study. 2003.

Smidt C.R. and Shieh D. Non-invasive biophotonic assessment of skin careotenoids as a biomarker of human antioxidant status. FASEB J 2003; submitted

Smidt, C.R. Ph.D., FACN, Effect of LifePak[®] Supplementation on Antioxidant Status Using BioPhotonic Raman Spectroscopy. Pharmanex in-house Study. 2002.

Stahl, W., Nicolai, S., Briviba, K., Hanusch, M., Broszeit, G., Peters, M., Martin, H.D. and Sies, H. (1997) Biological activities of natural and synthetic carotenoids: induction of gap junctional communication and singlet oxygen quenching. *Carcinogenesis* **18**, 89-92.

Steinberg FM, Chait A. Antioxidant vitamin supplementation and lipid peroxidation in smokers. *Am J Clin Nutr.* 1998 Aug;68(2):319-27.

Svilaas A, Sakhi AK, Andersen LF, Svilaas T, Strom EC, Jacobs DR Jr, Ose L, Blomhoff R. Intakes of antioxidants in coffee, wine, and vegetables are correlated with plasma carotenoids in humans. *J Nutr.* 2004 Mar;134(3):562-7.

Trosko, J.E. (2003) The role of stem cells and gap junctional intercellular communication in carcinogenesis. *Journal of Biochemistry and Molecular Biology.* **36**, 43-48.

Tsuchiya M., Scita, G., Thompson, D.F.T., Packer, L., Kagan V.E., Livrea, M.A. Retinoids and Carotenoids are peroxyl radical scavengers in Retinoids-progress in research and clinical applications, Livrea, M.A. and Packer, L (eds) Marcel Dekker Inc. New York, 1993.

Vogt, T.M., Mayne, S.T., Graubard, B.I., Swanson, C.A., Sowell, A.L., Schoenberg, J.B., Swanson, G.M., Greenberg, R.S., Hoover, R.N., Hayes, R.B. and Zeigler, R.G. (2002) Serum lycopene, other serum carotenoids, and risk of prostate cancer in US Blacks and Whites. *American Journal of Epidemiology* **155**, 1023-1032.

Yale Cancer Center News Releases: *Yale Researcher Funded to Study Nutritional Biomarkers.* May 19, 2003. http://www.yalecancercenter.org/ycc/releases/030519_35.html

Yamasaki, H. (1995) Non-genotoxic mechanisms of carcinogenesis: studies of cell transformation and gap junctional intercellular communication. *Toxicology Letters.* **77**, 55-61.

Yamasaki, H., Mesnil, M., Omori, Y., Mironov, N. and Krutovskikh, V. (1995) Intercellular communication and carcinogenesis. *Mutation Research* **333**, 181-188.

Zhang, L.X., Cooney, R.V. and Bertram, J.S. (1991) Carotenoids enhance gap junctional communication and inhibit lipid peroxidation in C3H/10T1/2 cells: relationship to their cancer chemopreventive action. *Carcinogenesis* **12**, 2109-2114.

Zhao, D.Y. et al, Resonance Raman measurement of macular carotenoids in retinal, choroidal, and macular dystrophies. *Arch Ophthalmol.* 2003 Jul;121(7):967-72.
http://www.ncbi.nlm.nih.gov/entrez/query.fcgi?cmd=Retrieve&db=pubmed&dopt=Abstract&list_uids=12860799